



# Purilex<sup>®</sup> Powder

(Cefalexin Monohydrate)



## Therapeutical category:

Oral broad-spectrum cephalosporin



## Characteristics:

Granules obtained by means of dry compaction of the crystalline powder without the addition of any excipients. No solvents are used in the process for higher purity with the exception of acetone used for the final product wash.



## Applications:

The product is suitable for the manufacture of granules for suspensions, dry powder mixtures for syrups and tablets, after the addition of appropriate excipients.

### This product is part of our PureActives<sup>®</sup> range.

All our PureActives<sup>®</sup> antibiotic products are produced using our proprietary, sustainable, and environmentally friendly enzymatic technology. They bring our brand promise to life through superior quality, outstanding reliability and leading sustainability performance.

All our manufacturing sites use the best available technology and operate dedicated waste water treatment plants, 24/7 throughout the year. Effluents are regularly sampled and checked for antimicrobial activity.



## Key parameters: Purilex<sup>®</sup> Powder (Cefalexin Monohydrate)

<b>Pharmacopeia quality</b>	USP, EP, BP	
<b>Regulatory information</b>	USA: FDA EU:EDQM Canada: HPB Japan: PMDA India: MoH South Korea: KFDA	Australia: TGA China: CFDA Japan: PMDA Taiwan: DoH Bangladesh: Moh Brazil: ANVISA
<b>Appearance</b>	White to almost white crystalline powder	
<b>Assay (anhydrous basis)</b>	Cephalexin: 97.0–102.0%	
<b>pH</b>	4.0 - 5.5	
<b>Water</b>	4.0 - 8.0%	
<b>Tapped bulk density</b>	≥ 0.65 g/ml	
<b>Total impurities</b>	≤ 3.0%	
<b>Residual solvents</b>	Acetone <0.2% (washing step)	
<b>Dimethylaniline</b>	Not used in the process	
<b>Shelf-life</b>	Minimum 5 years (prod. DSP Spain) or 4 years (prod. DSP China) in the original packaging under storage conditions	
<b>Batch size</b>	Approximately 1000 - 2400 kg (depending on manufacturing site)	
<b>Storage conditions</b>	Below 25°	
<b>Packaging</b>	Primary packaging: polyethylene bag	
	Secondary packaging: thermo-sealed aluminum polyethylene bag	
	Outer packaging: corrugated box or HDPE drum	
	Fixed quantity of 25 kg packed in each bag	

### Contact and information

For more information, please visit [www.dsm-sinochem.com](http://www.dsm-sinochem.com) or e-mail: [info@dsm-sinochem.com](mailto:info@dsm-sinochem.com) or contact us at one of the addresses below:

**Singapore (HQ)**  
sales@dsm-sinochem.com

**Europe/North America**  
Delft, The Netherlands  
sales.ea@dsm-sinochem.com

**Asia Pacific / Middle East / Africa**  
Gurgaon, India  
sales.amea@dsm-sinochem.com

**China**  
Beijing, P.R. China  
sales.china@dsm-sinochem.com

**Mexico & Latin America**  
Ramos Arizpe, Mexico  
sales.mla@dsm-sinochem.com

**Copyright © 2016 DSM Sinochem Pharmaceuticals.** All rights reserved. Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty as to the accuracy, currency or completeness of this information nor can it be construed to imply any other representation or warranty of any kind, including the implied warranties of merchantability, fitness for a particular purpose or non-infringement.

The content of this document is subject to change without further notice. Please contact us for the latest version of this document or further information. Product offers are explicitly not made in respect of those jurisdictions where valid third party patent rights related to such products are in force. In such jurisdictions the advertised products are only offered for regulatory purposes insofar allowed under local provisions like the Bolar provisions in the US and 2004/27/EC in Europe.  
VAPFS 2016